



**MedImmune**

AstraZeneca 

## **Developing Highly Productive Bioprocesses to Prepare for Pandemic Outbreaks**

**Jonathan Liu and Mark Thompson**

**International Vaccine Technology Workshop**

**Hyderabad, India**


**September 18, 2010**

# Influenza Vaccines

- Vaccination is the most effective way to prevent infection  
(CDC, 2009 <http://www.who.int/mediacentre/factsheets/fs211/en/index.html>)
- Two types of vaccines available
  - ◆ Inactivated vaccine – i.m.
  - ◆ Live attenuated influenza vaccine (LAIV) - Intranasal
- MedImmune vaccine (LAIV) approved for use in
  - ◆ USA
  - ◆ Canada
  - ◆ Korea
  - ◆ HK and
  - ◆ other countries




# Influenza Vaccines and Reverse Genetics (I)




This new technology eliminates the risk from potential wild type virus contaminants

# Influenza Vaccines and Reverse Genetics (II)





# Flu Vaccines and Cell Culture-based Production



|   | PRODUCTION SUBSTRATE            |              |
|---|---------------------------------|--------------|
|   | Eggs (SPF)                      | Cell Culture |
| Exposure of flock to environmental agents | Low risk, high impact           | NA           |
| Preproduction characterization            | - Limited                       | - Extensive  |
|   | - Inherent contamination        | - Sterile    |
| Manufacturing procedures                  | - Need to control contamination | Controlled   |
| Egg allergies limit use                   | Yes                             | No           |
| Process scalability                       | Slow and difficult              | Simple       |

**Cell culture-based vaccines provide a robust production platform in pandemic outbreak preparedness**

# MedImmune MDCK Cells - A Selected Production Cell Substrate



- A MDCK cell clone was selected from 13 cell lines based on**
- ✓ High productivity
  - ✓ Serum free growth
  - ✓ Susceptibility to a wide range of flu virus infections
  - ✓ Best cell clone out of > 2500 clones (refer to next slide)

# MedImmune MDCK Cells - Well Characterized Cell Banks

| <i>Test</i>  | <i>Cell Seed</i> | <i>MCB</i> | <i>WCB</i> | <i>ECB<sup>1</sup></i> |
|--|------------------|------------|------------|------------------------|
| <b>1. IDENTITY AND PURITY</b>  |                  |            |            |                        |
| Morphology   | +                | +          | +          | +                      |
| Identification e.g., isoenzymes, immunological and cytogenetic markers, DNA fingerprinting | +                | +          | +          | +                      |
| Karyotype  | +                | +          | -          | +                      |
| Life span (diploid cell lines)   | -                | -          | -          | +                      |
| Viability  | -                | +          | +          | +                      |
| Genetic stability (engineered cell lines)  | -                | +          | -          | +                      |
| <b>2. EXTRANEIOUS AGENTS</b>   |                  |            |            |                        |
| Sterility <sup>2</sup>   | +                | +          | +          | +/-                    |
| Electron microscopy  | -                | +          | -          | +                      |
| Tests in cell cultures   | -                | +          | +          | +                      |
| Retroviruses <sup>3</sup>  | -                | -          | -          | +                      |
| Tests in animals and eggs  | -                | +          | +          | +                      |
| Selected viruses by molecular methods  | -                | -          | -          | +                      |
| Antibody-production tests (rodent cell lines <sup>4</sup> )                                | -                | -          | -          | +                      |
| Bovine and/or porcine viruses  | -                | -          | -          | +                      |
| <b>3. BIOLOGICAL CHARACTERISTICS</b>   |                  |            |            |                        |
| Growth Characteristics   | -                | +          | +          | +                      |
| Tumorigenicity <sup>5</sup>  | -                | -          | -          | +                      |
| Oncogenicity <sup>6</sup>  | -                | -          | -          | +                      |

- ✓ No detectable adventitious agents
- ✓ No detectable tumorigenicity and oncogenicity
- ✓ Stable beyond End of Production level

WHO, 4 May 2010


# MedImmune MDCK Cells

## - Similarity of viruses with Commercial LAIV

| Analytical Test                                 | Comparability between egg and cell produced vaccine |
|---|---|
| Complete Genomic Sequence                       | ✓   |
| Phenotypic Analysis ( <i>ca</i> and <i>ts</i> ) | ✓   |
| Host Cell Susceptibility                        | ✓   |
| Virus Protein Expression                        | ✓   |
| Virus Morphology and Size                       | ✓   |
| Replication and Attenuation in Ferrets          | ✓   |
| Immunogenicity and Efficacy in Ferrets          | ✓   |
| Safety profile in Animal Models                 | ✓   |




# Disposable Culture Vessel-based Manufacturing Process



## ■ Fully disposable process implemented in GMP Pilot Plant

- ◆ No need for cleaning/validation with disposable culture vessels
- ◆ Shortened timeline for implementation in clinical production


## ■ Quick turned-around between batches (in a few hours) making possible to re-start production very rapidly



■ Also fully developed in MedImmune


# Annual Flu Vaccine Production Timeline

|                      | JAN | FEB | MAR | APR | MAY | JUN | JUL | AUG | SEP | OCT | NOV | DEC |
|----------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <b>Surveillance</b>  |     |     |     |     |     |     |     |     |     |     |     |     |
| <b>Selection</b>     |     |     |     |     |     |     |     |     |     |     |     |     |
| <b>Reassortment</b>  |     |     |     |     |     |     |     |     |     |     |     |     |
| <b>Manufacturing</b> |     |     |     |     |     |     |     |     |     |     |     |     |
| <b>Distribution</b>  |     |     |     |     |     |     |     |     |     |     |     |     |
| <b>Immunization</b>  |     |     |     |     |     |     |     |     |     |     |     |     |




- Frequent vaccine strains and manufacture process changes
- Short (often 1-3 weeks) Process Development Time

# Platform Production Reduces Production and Development Time



- Productivity increased 3 - 50 folds through process optimization
- Process robustness improved by identifying process variables and optimizing critical parameters

# Low Level of Impurities



Final Bulk DNA Level :  $\leq 1$  ng per dose

# Cell Culture vs. Egg-based Production of Flu Vaccine

| Production Time in Weeks* | MDCK Cells** | Eggs*** |
|---------------------------|--------------|---------|
| 4                         | 112          | 26      |
| 8                         | 223          | 52      |
| 12                        | 336          | 78      |

\* million doses of blended vaccine bulk produced

\*\* at 2 x 2000L scale excluding 3 wk lead time (cell thaw and seed train)

\*\*\* at 32,500 egg scale excluding lead time (egg procurement)

**Large amount of vaccine bulk can be produced in much shorter time using MDCK cells**

## ■ Intranasal administration

- ◆ simplifies the mass immunization process

## ■ MDCK cell–based production technology

- ◆ can produce flu vaccines quickly and in large quantity

## ■ Disposable technology

- ◆ minimizes change-over time between runs and leads to significant upfront savings

# Acknowledgements of US Government Sponsorship

This project is funded in whole or in part with Federal funds from the Office of the Assistant Secretary for Preparedness and Response (ASPR), Biomedical Advanced Research and Development Authority, under Contract Nos. HHSO100200600010C and HHSO100200700036C. The total federal program funding for these contracts is \$221,379,570, representing approximately 92% of the total amount of the projects. The remaining 8% of the total amount for the projects is anticipated to be financed by nongovernmental sources

The views expressed do not necessarily reflect the views or official policies of the Department of Health and Human Services; nor does mention of trade names, commercial practices, or organizations imply endorsement by the U.S. Government



# Questions



**Telephone No.**  
**1-650-603-2576**

**E-mail Address**  
**[liuj@medimmune.com](mailto:liuj@medimmune.com)**